
Lee's Cerclage Balloon (Uniconcave Balloon)

Lee's Cerclage Balloon is used for repositioning the prolapsed fetal membrane into the uterus prior to performing cervical cerclage to prevent still birth or preterm delivery

§ Features of Products

- Soft silicone balloon reduces risk of rupture of fetal membrane during procedure.
- Range of proper air inflation volume and diameter of inflated balloon printed on product
- Depth marking printed on product for easier measurement of position of balloon and fetal membrane

Introduced on cover of January 2015 issue of *American Journal of Obstetrics & Gynecology*

Prod. No	Air Fill Volume (Balloon should be completely deflated prior to inflation)	Diameter of Inflated Balloon	Diameter of Prolapsed Fetal Membrane
AI-001	3 ~ 4 mL	10 ~ 16 mm	Less than 10mm
AI-002	6 ~ 10 mL	17 ~ 23 mm	Over 10mm ~ Less than 20mm
AI-003	10 ~ 15 mL	25 ~ 30 mm	Over 20mm ~ Less than 25mm
AI-004	15 ~ 30 mL	30 ~ 40 mm	Over 25mm

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Direction for Use

1. Select suitable size of device.
2. Fill suggested air volume into balloon with needleless syringe, and check for air leakage from balloon surface or balloon bonding area (or check valve).
After leakage test, deflate balloon completely by releasing syringe plunger.
3. Procedure is performed under general anesthesia in Trendelenburg position.
4. Transabdominal amniocentesis is performed under ultrasonographic guidance, to reduce amount of amniotic fluid.
5. Fetal membrane and cervix are exposed using two Simpson retractors, and cervical edges are gently grasped with two atraumatic forceps.
6. Fill suggested air volume into balloon with needleless syringe, and insert device until balloon is properly positioned on prolapsed fetal membrane.
7. Gently push device until prolapsed fetal membrane is repositioned into uterus.
8. Perform cerclage procedure.
9. After cerclage procedure, balloon is deflated and device is removed from the patient.

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Yushin Medical & Trading Co., Ltd.