


Urology





Foley Balloon Catheters, 2 – Way, Silicone, Standard Type

Prod. No	Size	Balloon Capacity	Length(mm)
1001-06	6Fr	1.5ml	330
1001-08	8Fr	3ml	
1001-10	10Fr		
1002-12	12Fr	5ml	420
1002-14	14Fr		
1002-16	16Fr		
1002-18	18Fr		
1002-20	20Fr		
1002-22	22Fr		
1002-24	24Fr		
1002-26	26Fr		
1004-14	14Fr	30ml	420
1004-16	16Fr		
1004-18	18Fr		
1004-20	20Fr		
1004-22	22Fr		
1004-24	24Fr		
1004-26	26Fr		

Foley Balloon Catheters, 2 – Way, Silicone, Tiemann Type

Prod. No	Size	Balloon Capacity	Length(mm)
1011-10	10Fr	3ml	420
1011-12	12Fr	5ml	
1011-14	14Fr		
1011-16	16Fr		
1011-18	18Fr		
1011-20	20Fr		
1011-22	22Fr		
1011-24	24Fr		




Foley Balloon Catheters, 3 – Way, Silicone, Standard Type

Prod. No	Size	Balloon Capacity	Length(mm)
1005-14	14Fr	30ml	420
1005-16	16Fr		
1005-18	18Fr		
1005-20	20Fr		
1005-22	22Fr		
1005-24	24Fr		

UniBal[®] Foley Balloon Catheters, 2 – Way, Silicone, Advanced Integrated Balloon Type

Prod. No	Size	Balloon Capacity	Length(mm)
1002 UB-12	12Fr	5ml	420
1002 UB-14	14Fr		
1002 UB-16	16Fr		
1002 UB-18	18Fr		
1002 UB-20	20Fr		
1002 UB-22	22Fr		
1002 UB-24	24Fr		
1002 UB-26	26Fr		



(Conventional Integrated Balloon)

(UniBal[®])

- UniBal[®] has no bonding parts (A&B) on the border of balloon. Therefore distal part of the shaft are really flat and smooth.
- UniBal[®] is a registered trademark of **Yushin Medical Co., Ltd.**



Nephrostomy Catheters, Silicone, Standard Type

Prod. No	Size	Balloon Capacity	Length(mm)
1021-08	8Fr	1.5ml	420
1021-10	10Fr		
1021-12	12Fr	2ml	
1021-14	14Fr	3ml	
1021-16	16Fr	5ml	
1021-18	18Fr		
1021-20	20Fr		
1021-22	22Fr		
1021-24	24Fr		

Nephrostomy Catheters, Silicone, Integrated Balloon Type

Prod. No	Size	Balloon Capacity	Length(mm)
1021 F-08	8Fr	1.5ml	400
1021 F-10	10Fr		
1021 F-12	1021 UB-12	2ml	
1021 F-14	1021 UB-14	3ml	
1021 F-16	1021 UB-16	5ml	
1021 F-18	1021 UB-18		
1021 F-20	1021 UB-20		
1021 F-22	1021 UB-22		
1021 F-24	1021 UB-24		



Suprapubic Catheters, Silicone, Integrated Balloon Type

Prod. No	Size	Balloon Capacity	Length(mm)
1031-10	10Fr	3ml	400
1032-12	1032 UB-12	5ml	
1032-14	1032 UB-14		
1032-16	1032 UB-16		
1032-18	1032 UB-18		
1032-20	1032 UB-20		
1032-22	1032 UB-22	22Fr	


Suprapubic Catheter Sets (Catheter, Trocar, Sheath Introducer, Scalpel & Plug With Cap)

Prod. No	Size	Balloon Capacity	Length(mm)
1033-10	10Fr	3ml	400
1034-12	1034 UB-12	5ml	
1034-14	1034 UB-14		
1034-16	1034 UB-16		
1034-18	1034 UB-18		

Suprapubic Catheters, Silicone, Pigtail Tip, Integrated Balloon Type

Prod. No	Size	Balloon Capacity	Length(mm)
1040-12	12Fr	5ml	460
1040-14	14Fr		
1040-16	16Fr		
1040-18	18Fr	10ml	

Remark	
1021 F, 1031, 1032 1033, 1034, 1040	Conventional Integrated balloon type
1021 UB, 1032 UB 1034 UB	UniBal ® type



- Integrated balloon is basically applied for suprapubic catheters.



Hemostatic / Continuous Irrigation Catheters, Silicone

Dufour Pattern

Prod. No	Size	Balloon Capacity	Length(mm)
1037-18	18Fr	50-80ml	400
1037-20	20Fr		
1037-22	22Fr	80-100ml	
1037-24	24Fr		

Mercier Pattern

Prod. No	Size	Balloon Capacity	Length(mm)
1038-18	18Fr	50-80ml	400
1038-20	20Fr		
1038-22	22Fr	80-100ml	
1038-24	24Fr		

Couvelaire Pattern

Prod. No	Size	Balloon Capacity	Length(mm)
1039-18	18Fr	50-80ml	400
1039-20	20Fr		
1039-22	22Fr	80-100ml	
1039-24	24Fr		



Dufour

Mercier

Couvelaire

- Different balloon capacity will be available upon request.



Clean Intermittent Catheterization Sets, Silicone (C.I.C or Selfcath)

- Clean Intermittent Catheter is supplied in plastic case.
- Disinfectant prescribed by physician should be filled in plastic case.
- Available in male, female, pediatric and Tiemann type.
- Supplied with carrying case.

Prod. No	Type	Fr	I.D mm	O.D mm	Length mm
1035 F-10	Female	10	1.7	3.3	165
1035 F-12	Female	12	2.0	4.0	165
1035 F-14	Female	14	2.5	4.7	165
1035 M-10	Male	10	1.7	3.3	300
1035 M-12	Male	12	2.0	4.0	300
1035 M-14	Male	14	2.5	4.7	300
1035 PF-09	Pediatric, Female	9	1.5	3.0	150
1035 PM-09	Pediatric, Male	9	1.5	3.0	265
1035 ML-12	Male, Long	12	2.0	4.0	350
1035 T-12	Tiemann	12	2.0	4.0	300

- Additional sizes will be available upon request.